OSSIFY

Allograft Particulate

OSSIFY™ ALLOGRAFT PARTICULATE | Surgical Esthetics®

MINERALIZED CORTICAL-CANCELLOUS MIX

- Mixture mimics native composition of human bone
- Cancellous (~30–40%) provides ideal porosity for rapid clot formation
- Cortical (~60–70%) provides excellent space maintenance with slower resorption

- **0.6–1.25mm Particle Size:** GBR, extraction socket, implant placement, ridge split.

CODE	BG-012 025	BG-012 050	BG-012 100	BG-012 200
VOLUME	0.25cc	0.5cc	1.0cc	2.0cc

- **1–2 mm Particle Size:** GBR, large defect, crestal or lateral sinus elevation, ridge augmentation, adjacent extraction sockets.

CODE	BG-013 050	BG-013 100	BG-013 200
VOLUME	0.5cc	1.0cc	2.0cc

MINERALIZED CANCELLOUS

- Ideal inner surface area and porosity allows for immediate, stable clot formation
- Mineralized bone provides osteoconductive scaffold for efficient turnover to vital host bone
- Optimal particle structure and shape for natural and complete bone fill

- **0.25–1mm Particle Size:** Extraction socket, implant placement, perio defect, ridge split, apicoectomy

CODE	BG-010 050	BG-010 100	BG-010 200
VOLUME	0.5cc	1.0cc	2.0cc

- **1–2 mm Particle Size:** Large defect, sinus elevation, ridge augmentation, adjacent extraction sockets

CODE	BG-011 050	BG-011 100	BG-011 200
VOLUME	0.5cc	1.0cc	2.0cc

MINERALIZED CORTICAL

- Mineralized bone provides osteoconductive scaffold for efficient turnover to vital host bone
- Cortical bone provides excellent space maintenance with slower resorption
- Optimal particle density for long-term volume maintenance

- **0.25–1mm Particle Size:** Large defect, GBR, crestal or lateral sinus elevation, ridge augmentation

CODE	BG-014 050	BG-014 100	BG-014 200
VOLUME	0.5cc	1.0cc	2.0cc



PASSIONATE ABOUT REGENERATION


clinic@surgicalesthetics.com 833-701-5245 www.surgicalesthetics.com


Surgical Esthetics®

PREMIUM SUTURES

DEMINERALIZED CANCELLOUS SPONGE

- Demineralized bone promotes osteoinductive BMP activity for ectopic bone generation
- Compresses during placement, yet expands to maintain space
- Rapidly absorbs blood for quick, stable clot formation
- Unlike block grafts, Sponge Strip won't cause exposures or soft tissue perforation





- **Sponge Strip:** GBR, vertical or horizontal ridge augmentation, extraction socket, implant placement, perio defect

CODE	BG-024010	BG-024020
SIZE	10x6x7mm	20x6x7mm

- **1-4mm Particulate:** Extraction socket, implant placement, patients with compromised bleeding

CODE	BG-025 075	BG-025 100
VOLUME	0.75cc	1.0cc



DEMINERALIZED CORTICAL

- Demineralized bone can promote osteoinductive BMP activity
- Excellent for small periodontal defects
- Optimal particle density for clot stability

- **0.25-1mm Particle Size:** Large defect, GBR, crestal or lateral sinus elevation, ridge augmentation

CODE	BG-020 050	BG-020 100	BG-020 200
VOLUME	0.5cc	1.0cc	2.0cc



PASSIONATE ABOUT REGENERATION

clinic@surgicalesthetics.com 833-701-5245 www.surgicalesthetics.com